


Poc-KIT stretcher diagram




The PocKIT Stretcher can be deployed in seconds and saves vital minutes in emergency casualty evacuation. To use the Poc-KIT Stretcher just follow the instructions below making sure to adhere with the principles of First Aid before moving the person.

PLEASE NOTE: The POC-KIT Stretcher is designed for emergency evacuation and must be used in accordance with the principles of First Aid.











There are grab handles points for up to 3 people to carry the stretcher when needed. Weight tested to carry up to 150kg (23 stone).

The Poc-KIT stretcher can also be used as a Bivvy. To set up the Bivvy just follow instructions below.

